


IMMEDIATE GRAFT


- Ready to use
- No mixing required



STORAGE

REFRIGERATE PRODUCT UPON ARRIVAL

Products are shipped in an insulated shipping box with cold packs. The cold packs may become thawed during shipment. Immediate Graft contains an organic molecule that stimulates bone formation that is temperature sensitive but stable at room temperature. If the small temperature indicator located on the product packaging is black upon arrival, the efficacy of the graft material may have been compromised. Notify us immediately if the dot is black. Refrigerate the product if the dot has not turned black.



DIRECTIONS

- After extraction and placement of immediate implant, fill voids between the implant and bone with Immediate Graft.
- After placing Immediate Graft, cover with a d-PTFE (Teflon) membrane or temporary crown. When using crowns, place the temporary abutment first grafting. After grafting is complete, cement temporary crown. **Do not use animal or human membranes to cover Immediate Graft or the graft may fail.**
- Place sutures as needed. Oral Bond is advised to be placed over the membrane, creating a hermetic seal to prevent oral fluids and bacteria from entering the grafted site.
- Systemic antibiotics are advised.

To avoid displacing the graft material, instruct patient to avoid the area with their tongue, and when eating or brushing. No vigorous rinses or spitting.

PRECAUTIONS

Do not leave defect open. Do not compromise blood supply to the defect area. This product has not been studied in pediatric patients.

Federal law restricts this device to sale by or on the order of a dentist.



See reverse for additional instructions

IMMEDIATE GRAFT

DESCRIPTION

Immediate Graft is a combination of Socket Graft and OsseoConduct β TCP granules (250-500 μ m). The β TCP granules are premixed in this product. This product contains no human or animal tissue.

SOCKET GRAFT

Socket Graft is a dual phase calcium phosphate based synthetic bone graft material. It is designed to quickly and safely treat the extraction socket, to retain alveolar bone, stimulate bone formation, and reduce post-operative complications. This product is designed to be quickly and completely resorbed, to produce normal, healthy vital bone.

Socket regeneration (preservation) insurance code: D7953

INDICATIONS

Immediate Graft is a bone graft material indicated for grafting small bony defects.

CONTRAINDICATIONS

This product is contraindicated in individuals who are unable to follow post-op instructions.

OSSEOCONDUCT

OsseoConduct is a synthetic, microporous, biocompatible, third generation beta tricalcium phosphate bone grafting material. OsseoConduct is resorbed and replaced by bone in 4-6 months and is intended for use as a general bone graft material. OsseoConduct is designed to maintain shape and volume when used for ridge augmentation or socket grafting and lateral window sinus augmentations. OsseoConduct is not approved for immediate support of removable dental appliances.

INDICATIONS

OsseoConduct is intended for use as a maxillofacial bone graft material in bone voids and for appositional growth.

CONTRAINDICATIONS

OsseoConduct is not indicated where loading of the graft site is planned, such as oral sites that are intended to be loaded immediately with removable dental prosthesis.

WARNINGS

Immediate Graft is single use only. Do not re-sterilize or reuse. Do not use Immediate Graft if the package has been opened or damaged, or has exceeded the expiration date. Do not use if the temperature indicator on the packaging has turned black.

STEINERBIO

1051 Olsen St., Bldg 3611
Henderson, NV 89011
contactus@steinerbio.com
www.steinerbio.com
866.317.1348